

ASSAY OF β GALACTOSIDASE IN INDUCED BACTERIAL CELLS page 35

March 1994, rvsd 27 Dec 1994, 3 Jan '97, 5 Jan '98, 7 Mar 00, 5 March 01, 3 Jan 03, 30Dec04, 31Dec09, 8Mar11

David Fankhauser, Professor of Biology and Chemistry

http://biology.clc.uc.edu/fankhauser/Labs/Genetics/Lac_Operon_Derepression/Lactase_assay_induced.htm

Your team should now have five samples of toluenized cells: 0, 20, 40, 60 and 80 minute lactose-grown from the previous stage of the experiment. Last week's preliminary experiment demonstrated the approximate repressed and derepressed levels of β galactosidase and optimum assay conditions for these toluenized cells. You will now assay the amount of β galactosidase in this set of cells using a similar protocol, and then calculate the specific activity to reveal derepression.

EQUIPMENT per desk of two:

15x 13 x 100 mm test tubes in TT rack
 1x 37°C hot block, thermometer
 1x 5 mL displacement pipetes & tips
 2x 1 mL micropipets & tips
 1x stopwatch
 1x spectrophotometer
 2x cuvettes, in TT rack

SUPPLIES (per 16 students):

your 5 samples of toluenized cells on ice:
 T-0, T-20, T-40, T-60, T-80
 (Minutes grown in lactose)
 350 mL 2% K₂CO₃ in 2.0 mL repipet
 20 mL 20 mM orthonitrophenyl galactoside
 50 mL 0.1 M PO₄ buffer, pH 6.0, 0.4 repeater
 25 mL dH₂O (in 125 mL flask)

1. Label tubes, per team of two:

five Rxn tubes: R-0, R-20, R-40, R-60, R-80
 ten K₂CO₃ tubes, start and finish: S0, S20, S40, S60, S80, F0, F20, F40, F60, F80

2. Add 2 mL 2% K₂CO₃ to each of the S and F tubes with repeater pipet.

3. Set up the assay table. Set up the following table in your notebook (2x or 4x spaced)

tube	sample	A ₆₆₀	dH ₂ O	buffer	toluenized cells	A ₄₁₅	F - S	$\Delta A_{415}/A_{660}/\text{aliquot}$
R-0	T-0	-----	--	0.4	3.2	S	_____	_____
						F	_____	_____
R-20	T-20	-----	- -	0.4	3.2	S	_____	_____
						F	_____	_____
R-40	T-40	-----	2.2	0.4	1.0	S	_____	_____
						F	_____	_____
R-60	T-60	-----	2.7	0.4	0.5	S	_____	_____
						F	_____	_____
R-80	T-80	-----	2.7	0.4	0.5	S	_____	_____
						F	_____	_____

4. **Set up reaction tubes:** add water, then buffer (with repeater pipet), then toluenized cells. Mix, place in 37°C hot block to pre-warm. Place a thermometer in tube R-0, wait until temp = 35°C. Because the large sample of cold cells in tubes R-0 and R-20, we will start the reaction in R-80 first, then R-60, etc. while tubes R-0 and R-20 continue to warm up.

5. Add 0.4 mL ONPG down the side of all five R tubes, vortex.

Then begin taking Start samples: at 30 second intervals:

- a: **Withdraw 1.0 mL start sample** from R-80 assay tube, add to the S-80 tube (contains K₂CO₃), mix, and **start stopwatch**. Place S-80 in rack to be read later. Wait 30 seconds.
- c: **Withdraw 1.0 mL R-60 tube**, add to S-60 exactly 30 seconds after the previous sample. Repeat for remaining reaction tubes, R-40, R-20 and R-0.
- d: **Incubate** the five assay tubes for exactly 15 minutes each at 37°C.
- e: **When the stopwatch reads 15 minutes, withdraw 1.0 mL from R-80**, add to F-80 (with K₂CO₃), mix. At 30 second intervals, repeat aliquoting to F-60, F-40, F-20 and F-0

4. **Read the A₄₁₅ of all tubes, reading S-0 to S-80 tubes first, then the F-0 to F-80 last.** Subtract the start from the finish = ΔA_{415} (which is proportional to the amount of enzyme present).

5. Calculate specific activities for each sample:

(1/mL of cells tested) x (ΔA_{415} from assay)/(A₆₆₀ of the culture at sampling time)

6. **Graph the specific activities of the two cultures versus time.** Discuss differences in specific activity seen as the culture grows for progressively longer times in lactose.